



Your support in life

Knee Braces Knee Supports Foot Care Products Knee
Bed Neck Collar Wrist Splints Truss For a Hernia Visc
chosis Arm&Supplies Corsets Body Orthosis Hip Ortho
w Supports Shoulder Supports Ankle Braces Visc

aurafix[®]
ORTHOPEDIC PRODUCTS

■ About Us

Our company, which was founded in 2002 and has been an expert in producing orthopedic ready-made products, orthoses, prostheses and rehabilitation products, established TSD Sağlık Ürünleri in the year 2007 in order to convey its accumulated knowledge into wholesale and marketing and began manufacturing orthopedic products in its own brand, aurafix, and became popular in a short period of time.

TSD Sağlık Ürünleri was founded upon a small workshop, then grew and developed thanks to clean labor, exclusive service and investments made for people in technology. Today, the company continues its activities with its great machinery. Promptly, it approaches its goal which is to be a global brand known both in Turkey and the world.

Including the production staff and R&D, all departments are in constant endeavor towards this goal, by following up with the latest innovations in the sector and producing various solutions for special needs.

We came to where we are thanks to you, by thinking about you as much as you do and we are dedicated to work for you without compromising our quality.

Your support in life...Aurafix....

Index

Neck Collars	4 - 6
Shoulder Supports	8 - 10
Arm Supports	11 - 13
Elbow Supports	14 - 17
Wrist Splints	18 - 23
Wrist Orthosis	24 - 25
Finger Splints	26 - 26
Corsets	28 - 38
Body Orthosis	39 - 41
Hernia Belts	42 - 42
Knee Supports	44 - 51
Knee Braces	52 - 54
Calf Protection and Tights	55 - 56
Ankle Supports	57 - 62
Ankle Orthosis	63 - 64
Hip Abduction Orthosis	65 - 67
Products in Standard Size	68 - 69
Elastic Supports	70 - 70
Foot Health Products	72 - 78
Medical Compression stockings and Arm Sleeve	80 - 86
Silicone Breast Form	88 - 88
Rehabilitation	90 - 90
Pillows and seat cushions	92 - 94
Stands & our certificates	95 - 96



Neck Collars



Reinforced Cervical Collar

Ref : 920

This is the sponge neck collar with special hard polyethylene support designed anatomically for the neck. The outer cover can be removed and washed.

Indications;

- 1-Cervical disc hernia,
- 2-Muscle injury,
- 3-Cervical arthritis,
- 4-Muscle spasm.

Size;

S : 28 - 33 cm
M : 34 - 39 cm
L : 40 - 45 cm



Soft Cervical Collar

Ref : 910


This is the soft sponge neck collar, anatomically designed for the neck. By reminding the patient that he needs to refrain from excessive cervical activity, it provides sensory feedback and provides him warmth around the neck. The outer cover can be removed and washed.

Indications;

- 1-Cervical disc hernia,
- 2-Muscle injury,
- 3-Cervical arthritis,
- 4-Muscle spasm.

Size;

S : 28 - 33 cm
M : 34 - 39 cm
L : 40 - 45 cm



Cervical Collar without Chin Support

Ref: 925

This is the neck collar designed anatomically without the plastic jaw support. Height adjustment is possible with velcro fasteners.

Indications;

- 1-Cervical disc hernia,
- 2-Cervical arthritis,
- 3-Post-trauma (for protection).

Size;

XS (0) : 28 - 31 cm
S (1) : 32 - 35 cm
M (2) : 36 - 39 cm
L (3) : 40 - 43 cm
XL (4) : 43 - 46 cm
XXL(5) : 46 - 50 cm



Cervical Collar with Chin Support

Ref: 926

This is the plastic neck collar, anatomically designed for the neck. This collar is a more restrictive type with its jaw support, limiting neck movement. Height adjustment is possible with velcro fasteners.

Indications;

- 1-Cervical disc hernia,
- 2-Injury,
- 3-Post-operation and –trauma.

Size;

XS (0) : 28 - 31 cm
S (1) : 32 - 35 cm
M (2) : 36 - 39 cm
L (3) : 40 - 43 cm
XL (4) : 43 - 46 cm
XXL(5) : 46 - 50 cm





Wire-Frame Collar

Ref: 950

Designed anatomically for the shape of the neck, this collar is semi-rigid with the capability of flexion and extension adjustment.




Indications;

- 1-Cervical disc hernia,
- 2-Cervical arthritis,
- 3-Post-operation (for protection).

Size;

- M (2) : 33 - 39 cm
L (3) : 40 - 46 cm



Nelson Collar

Ref: 930

Manufactured with 12 mm perforated plastizode material.

Indications;

- 1-Cervical disc hernia,
- 2-Injury,
- 3-Post-trauma (for protection).

Size;

- XS (0) : 28 - 31 cm
S (1) : 32 - 35 cm
M (2) : 36 - 39 cm
L (3) : 40 - 43 cm
XL (4) : 43 - 46 cm



Nelson Collar with chin Support

Ref: 935

This has the same specifications as Ref 930. Differently, this collar is a more restrictive type with its jaw support, limiting neck movement.

Indications;

- 1-Cervical disc hernia,
- 2-Injury,
- 3-Post-trauma (for protection).

Size;

- XS (0) : 28 - 31 cm
S (1) : 32 - 35 cm
M (2) : 36 - 39 cm
L (3) : 40 - 43 cm
XL (4) : 43 - 46 cm





Philadelphia Collar

Ref: 940

This is the semi-rigid neck collar with front and rear support, manufactured with 12 mm perforated plastizode material.



Indications;

- 1-Soft tissue injury,
- 2-Stable bone and ligament injury,
- 3-Post-operation and -trauma (for protection and stabilization).

Size;

- XS : 28 - 32 cm
S : 33 - 37 cm
M : 38 - 42 cm
L : 42 - 46 cm



Philadelphia Tracheostomy Collar

Ref: 945

Used when Ref 940 can be used and additionally, in cases of tracheostomy.



Size;

- XS (0) : 28 - 32 cm
S (1) : 33 - 37 cm
M (2) : 38 - 42 cm
L (3) : 42 - 46 cm



Emergency Neck Collar

Ref: 948

This is the emergency collar manufactured with plastic material, with the inner side laminated.

Size;

- XS (0) : 28 - 32 cm
S (1) : 33 - 37 cm
M (2) : 38 - 42 cm
L (3) : 42 - 46 cm






Upper extremity



Clavicle Bandage

Ref: E-140

Manufactured with a sponge fabric with terry cloth inside. With the non-gradual velcro adjustment, it enables the shoulders to be pulled back as desired for a proper posture. Easy to use with its soft and sweat-preventing material.



Indications;

- 1-Clavicular fracture,
- 2-Clavicular and sternoclavicular dislocation.

Size;

XS : 18 - 24 cm
S : 25 - 30 cm
M : 31 - 36 cm
L : 37 - 42 cm
XL : 43 - 48 cm
XXL: 49 - 54 cm



Clavicle Bandage (with sponge)

Ref: E-142

Manufactured with cotton fabric over sponge. With the laces on the rear side, the shoulders can be pulled back as much as desired.



Indications;

- 1-Clavicular fracture,
- 2-Clavicular and sternoclavicular dislocation.

Size;

XS : 92 - 100 cm
S : 100 - 109 cm
M : 109 - 118 cm
L : 118 - 127 cm
XL : 127 - 136 cm
XXL: 136 - 145 cm





Acromio-clavicular Bandage

Ref: E-143

Manufactured with respirable perforated fabric. Shoulder and axillary pad is made of sponge-terry cloth for maximum comfort. The shoulder pad allows the load to be distributed evenly. The rotating buckle system enables easy wear.



Indications;

1-Acromioclavicular
dislocation- subluxation.

Size;

Left / Right Available.

S : 30 - 35 cm

M : 35 - 40 cm

L : 40 - 45 cm



Shoulder Support

Ref: 700

Made of neoprene material. The compression it applies prevents edema and reduces pain by decreasing the muscle tonus. The velcro fasteners at the bottom end provide ease of use. There is an axillary pad.

Indications;

- 1-Subluxation,
- 2-Tendinitis, osteoarthritis,
- 3-After getting a shoulder endoprosthesis,
- 4-Post-operation.

Size;


S : 22 - 24 cm

M : 24 - 27 cm

L : 27 - 30 cm

XL : 30 - 34 cm

XXL: 34 - 38 cm





Humerus Brace

Ref: 715

This is a shoulder orthosis made of rigid sponge. It is used for shoulder dislocations, and also for protecting and stabilizing the shoulder after operation and encasement in plaster. Use and application is quite simple. With the optional arm sling, the hand does not hang and edema does not form in the distal.



Indications;

- 1-Stable proximal fractures of the head of the humerus,
- 2-Post-operative and encasement in plaster,
- 3-For protection and stabilization in shoulder dislocation.

Size;

S : 22 - 24 cm
M : 24 - 27 cm
L : 27 - 30 cm
XL : 30 - 34 cm
XXL: 34 - 38 cm



Woolen Shoulder Support (Double Sided)


Ref: 705

Made of wool, it provides warmth and compression.



Size;

S : 22 - 24 cm
M : 24 - 27 cm
L : 27 - 30 cm
XL : 30 - 34 cm
XXL: 34 - 38 cm





Velpau Bandage

Ref: AO-01

Made of sponged fabric. Provides immobilization of the shoulder and the upper arm in the adduction. The bandages which go around the chest increase stability.



Indications;

- 1-Shoulder dislocation,
- 2-Post-op treatment of proximal fractures of the head of the humerus and scapula.
- 3-Rotator cuff lesions.

Size;

S : 30 - 33 cm
M : 33 - 36 cm
L : 36 - 39 cm
XL : 39 - 42 cm
XXL: 42 - 45 cm



Arm Sling

Ref: AO-02

Manufactured with a sponge fabric with terry cloth inside. The waist-supported strap provides extra protection.



Indications;

- 1-Shoulder dislocation-subluxation,
- 2-Post-operative for immobilization,
- 3-Prevention of edema forming in the distal in case of fracture.

Size;

S : 30 - 33 cm
M : 33 - 36 cm
L : 36 - 39 cm
XL : 39 - 42 cm
XXL: 42 - 45 cm





Arm Sling (Standard)

Ref: AO-02K

This is an arm sling used for immobilizing the shoulder and the upper arm, made of sponged fabric. Comes in standard size.

Size;

Standard Size.



Arm Sling (Simple)

Ref: AO-02BZ

This is a simple arm sling ,made of sponge fabric. Comes in standard sizes.

Size;


Standard Size.



Arm Sling (Perforated)

Ref: AO-02H

Manufactured with respirable and sweat-preventing perforated fabric. Used for the same cases as Ref AO-02.



Size;

S : 30 - 33 cm
M : 33 - 36 cm
L : 36 - 39 cm
XL : 39 - 42 cm
XXL: 42 - 45 cm





Arm Sling (30 Degrees)

Ref: AO-03

Manufactured with a sponge fabric with terry cloth inside. AO-03 enables a 30-degree immobilization in abduction. Prevents post-operational cohesion. Helps the patient in mobilization and rehabilitation. Relieves subacromial pressure.

Indications;

- 1-Post-operative impingement syndrome,
- 2-After rotator cuff rupture repair,
- 3-After manipulation of the frozen shoulder under anesthesia,
- 4-After shoulder endoprosthesis.

Size;

Standard Size.



Arm Sling (45-60 Degrees)

Ref: AO-04

Manufactured with a sponge fabric with terry cloth inside. AO-04 optionally enables a 45- or 60-degree immobilization in abduction. Prevents post-operational cohesion. Helps the patient in mobilization and rehabilitation. Relieves subacromial pressure. The hand exercise ball within the AO-04 arm sling is effective in accelerating the blood flow.



Indications;

- 1-Post-operative impingement syndrome,
- 2-After rotator cuff rupture repair,
- 3-After manipulation of the frozen shoulder under anesthesia,
- 4-After shoulder endoprosthesis.

Size;

Standard Size.





Tennis Elbow Support

Ref: 301

Made of neoprene material. It has a micro-massaging pad which has a droplet feature designed anatomically, which applies compression on the muscle groups adherent to the lateral epicondylitis.



Indications;

- 1-Lateral epicondylitis,
- 2-Increase of intra articular liquid,
- 3-Elbow pain.

Size;

Left / Right Available.

- XS : 20 - 22 cm
S : 22 - 24 cm
M : 24 - 26 cm
L : 26 - 28 cm
XL : 28 - 30 cm
XXL: 30 - 32 cm



Elbow Support with pad

Ref: 302

Manufactured anatomically with neoprene material. Because the perimeter of the elbow is supported with a laminated pad, it provides better protection.

Indications;

- 1-Protection against trauma,
- 2-Olecranon bursitis,
- 3-Protection during and after sports activities.

Size;

- XS : 20 - 22 cm
S : 22 - 24 cm
M : 24 - 26 cm
L : 26 - 28 cm
XL : 28 - 30 cm
XXL: 30 - 32 cm



Elastic Elbow Support

Ref: 300

Manufactured anatomically with neoprene material. The velcro fastener on the upper arm provides ease of use.

Indications;

- 1-Protection against trauma,
- 2-For providing heat and compression,
- 3-Protection during and after sports activities.

Size;

- XS : 20 - 22 cm
S : 22 - 24 cm
M : 24 - 26 cm
L : 26 - 28 cm
XL : 28 - 30 cm
XXL: 30 - 32 cm





Epicondylitis Bandage (with Gel Pad)

Ref: 309

This is an epicondylitis band with 3-point silicone pad for local compression and gel pad for cold treatment. The elastic band enables dynamic compression. The silicone pad can be switched among right and left elbows. Used for medial or lateral epicondylitis.



Size;
Standard Size.



Epicondylitis Bandage

Ref: 308

This is an epicondylitis band which has a 3-point silicone pad providing local compression. The elastic band enables dynamic compression. The silicone pad can be switched among right and left elbows. Used for medial or lateral epicondylitis.



Size;
Standard Size.



Neoprene Epicondylitis Bandage

Ref: 303

Manufactured anatomically with neoprene material. The 3-point pad providing local compression also reduces the tension in the tendon. The pad can be switched among right and left elbows. Used for medial or lateral epicondylitis.



Size;
Standard Size.





Woven Epicondylitis Elbow Support

Ref: 304

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. On both sides of the elbow, it has a micro-massaging silicone pad which has a droplet feature designed anatomically, which applies compression on the muscle groups adherent to the condyles. Local compression reduces tendon tension and thus pain.



Indications;

- 1-Medial, lateral epikondilitte,
- 2-Eklem içi sıvı artışında,
- 3-Ameliyat veya yaralanma sonrası iritasyonda kullanılır.

Size;

- S : 22 - 24 cm
M : 24 - 26 cm
L : 26 - 28 cm
XL : 28 - 30 cm
XXL: 30 - 32 cm



Woven Elbow Support

Ref: 306

Manufactured with respirable, flexible woven fabric, anatomically convenient for the shape of the elbow.

Indications;

- 1-Arthritis,
- 2-Increase of intra articular liquid,
- 3-Post-operation or -injury.

Size;

- S : 22 - 24 cm
M : 24 - 26 cm
L : 26 - 28 cm
XL : 28 - 30 cm
XXL: 30 - 32 cm



Epicondylitis Splint

Ref: AO-22

It is manufactured with respirable cotton fabric. It has an aluminum bar which restricts wrist motion. The fingers are free after the MCP joint.

Indications;

- 1-Epicondylitis,
- 2-Carpal tunnel syndrome,
- 3-In cases where the immobilization of the wrist is needed.

Size;

- S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm
XXL: 22 - 24 cm





Elbow R.O.M . Orthosis

Ref: 305

This is an adjustable elbow angle orthosis manufactured anatomically in accordance with elbow shape. Its made of sponge –fabric whose interior side is terry cloth. Angle adjustment can be fulfilled easily by a pressing buton. The orthosis provides controlled elbow flexion and extension. Angle flexion can be adjusted in degrees, (0-15-30-45-60-75-90-105) and extension in degrees (-15-0-15-30-45-60-75-90). According to the length of the arm, a length adjustment on the orthosis can be performed. There is a length-adjustable aluminum support preventing the wrist from falling. The shoulder strap provides ease of use.

Indications;

- 1-Helping rehabilitation in elbow contracture,
- 2-Providing protection and helping rehabilitation after fracture,
- 3-After ligament injury,
- 4-When limited and controlled motion of the elbow is desired,
- 5-Conservative treatment after elbow dislocation,
- 6-Providing protection and rehabilitation after elbow ligament instabilities in post operation period
- 7-Providing rehabilitation after elbow endoprosthesis in post operation period

Size;


Left / Right Available.
S : 30 - 37 cm
M : 38 - 45 cm



Cubital Tunnel Elbow Support

Ref:307

Made of neoprene material. Inside the elbow, there is a 20-degree inclined steel support. It prevents the elbow joint from bending thus ulnar nerve from straining and thus symptoms are relieved.



Indications;

- 1-Cubital tunnel syndrome,
- 2-Ulnar nerve problems.

Size;

XS : 20 - 22 cm
S : 22 - 24 cm
M : 24 - 26 cm
L : 26 - 28 cm
XL : 28 - 30 cm
XXL: 30 - 32 cm





Neoprene De Quervain Splint

Ref: 601


Made of neoprene material. It has an aluminum support restricting the thumb motion whose inclination can be adjusted.

Indications;

- 1-De Quervain tenosynovitis,
- 2-1.CMC and MCP joint irritation,
- 3-Protection in thumb ligament laxity,

Size;

Left / Right Available.
S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm




Bandage Wristband

Ref: 602

Made of neoprene material. It wraps around the wrist, providing a better protection. Used for protection and retaining warmth in wrist injuries, sprains and trauma.

Size;

Left / Right Available.
Standard Size.



Simple Wristband

Ref: 600

Simple wristband made of neoprene material. Used for protection and retaining warmth in wrist injuries, sprains and trauma.

Size;

Standard Size.





Neoprene De Quervain Splint with Thermoplastic Support Ref : 605

Thumb splint made of neoprene material. The support made of thermoplastic material which can be formed with hot water enables the CMC and MCP joint to be immobilized in the desired position.



Indications;

- 1-De Quervain tenosynovitis,
- 2-1.CMC and MCP joint degeneration,
- 3-After injury.

Size;

Left / Right Available.
Standard Size.



Neoprene Standard Wrist Splint Ref: 606

Wrist splint made of neoprene material. It has an anatomically designed support which provides the immobilization of the wrist in neutral position. Because the aluminum support is portable, the product can be used for the opposite hand. Comes in 2 sizes.



Indications;

- 1-Carpal tunnel syndrome,
- 2-Tendinitis, tenosynovitis, tenovaginitis,
- 3-For recovery in case of wrist pain,
- 4-For protection after injury and operation.

Size;

S : 14 - 18 cm
M : 18 - 22 cm





Neoprene Wrist Splint

Ref: 603

Wrist splint made of neoprene material. It has an anatomically designed support which provides the immobilization of the wrist in neutral position. There is a bandage around the wrist which increases stabilization.

Indications;

- 1-Carpal tunnel syndrome,
- 2-Tendinitis, tenosynovitis, tenovaginitis,
- 3-For immobilization in case of wrist pain,
- 4-For protection after injury and operation,
- 5-In fractures after encasement in plaster.

Size;

Left / Right Available.
S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm
XXL: 22 - 24 cm



Neoprene Wrist Splint with Thumb Support

Ref : 604

Made of neoprene material. It has an anatomically designed support which provides the immobilization of the wrist and thumb in neutral position. There is a bandage around the wrist which increases stabilization.




Indications;

- 1-In cases where Carpal tunnel and De Quervain are seen together,
- 2-De Quervain tenosynovitis,
- 3-1.CMC and MCP joint degeneration,
- 4-For immobilization in case of wrist and thumb pain,
- 5-For protection after injury and operation,
- 6-In fractures after encasement in plaster.

Size;

Left / Right Available.
S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm
XXL: 22 - 24 cm






Lace up Wrist Splint

Ref: 609

Manufactured with respirable perforated fabric. On the volar side of the wrist, there is an aluminum support which is anatomically formed. The bandage provides stabilization for the lace system around the wrist. Easy to use and apply. Used for the same cases as Ref-603. Comes in 3 sizes.



Size;

Left / Right Available.

S : 14 - 17 cm

M : 18 - 21 cm


L : 22 - 25 cm



Lace up Wrist Splint with Thumb Support

Ref: 610

Manufactured with respirable perforated fabric. It has an anatomically designed support which provides the immobilization of the wrist and thumb in neutral position. The bandage provides stabilization for the lace system around the wrist. Easy to use and apply. Used for the same cases as Ref-604. Comes in 3 sizes.




Size;

Left / Right Available.

S : 14 - 17 cm

M : 18 - 21 cm

L : 22 - 25 cm





Wrist Splint

Ref: AO-19

It is manufactured with respirable cotton fabric. On the volar side of the wrist, there is an aluminum support which is anatomically formed and provides the immobilization of the wrist in neutral position. The fingers are free after the MCP joint. There is a stabilization bandage around the wrist.


Indications;

- 1-Carpal tunnel syndrome,
- 2-Tendinitis, tenosynovitis, tenovaginitis,
- 3-For recovery in case of wrist pain,
- 4-For protection after injury and operation.

Size;

Left / Right Available.

- S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm
XXL: 22 - 24 cm



De Quervain Splint

Ref : AO-15

It is manufactured with respirable cotton fabric. It has an aluminum support restricting the thumb motion whose inclination can be adjusted.

Indications;

- 1-De Quervain tenosynovitis,
- 2-1.CMC and MCP joint irritation,
- 3-Protection in thumb ligament laxity.

Size;

Left / Right Available.

- S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm



Wrist Splint with Thumb Support

Ref : AO-21

It is manufactured with respirable cotton fabric. It has an anatomically designed support which provides the immobilization of the wrist and thumb in neutral position. There is a bandage around the wrist which increases stabilization.


Indications;

- 1-In cases where Carpal tunnel and De Quervain are seen together,
- 2-De Quervain tenosynovitis,
- 3-1.CMC and MCP joint degeneration,
- 4-For recovery in case of wrist and thumb pain,
- 5-For protection after injury and operation.

Size;

Left / Right Available.

- S : 14 - 16 cm
M : 16 - 18 cm
L : 18 - 20 cm
XL : 20 - 22 cm
XXL: 22 - 24 cm




Wrist Splints



Woven Wrist Splint

Ref: AO-17

Manufactured with flexible and sweat-preventing woven fabric. The silicone support around the ulnar styloid both protects the bone protuberance and reduces the pressure on veins and nerves. The volar side has a support providing immobilization of the wrist in neutral position. There is a stabilization bandage around the wrist.



Size;

Left / Right Available.

S : 14 - 16 cm

M : 16 - 18 cm

L : 18 - 20 cm

XL : 20 - 22 cm

XXL: 22 - 24 cm



De Quervain Splint

Ref: 611


This is a thumb splint which has a 1 mm polyethylene support around the thumb, stabilizing the 1 MCP joint. Used in 1 MCP joint irritation.

Size;

Left / Right Available.

M : 14 - 19 cm

L : 19 - 24 cm





Thermoplastic Hand Immobilization Splint

Ref : ORT-06

This is a hand and wrist orthosis made of thermoplastic material, covered with terry cloth inside. It keeps the wrist and the fingers in recovery position.



Indications;


- 1-After encasement of fractures,
- 2-For treatment of sprains and injuries,
- 3-Hemiplegia,
- 4-Post-operation,
- 5-For protection in spasticity deformity.

Size;

Left / Right Available.

M : 14 - 18 cm

L : 18 - 22 cm



Thermoplastic Wrist Immobilization Splint

Ref : ORT-05

This is a wrist orthosis made of thermoplastic material, covered with terry cloth inside. Keeps the wrist in neutral position. The fingers are free after the MCP joint.

Indications;


- 1-After encasement of fractures,
- 2-For treatment of sprains and injuries,
- 3-Carpal tunnel syndrome.

Size;

Left / Right Available.

M : 14 - 18 cm

L : 18 - 22 cm



Thermoplastic Anti-Spasticity Hand Splint

Ref : ORT-08


This is a hand and wrist orthosis made of thermoplastic material, covered with plastizode inside. Used especially for hand spasticity because of the grooved finger design and resistance to tension.

Size;

Left / Right Available.

M : 14 - 18 cm

L : 18 - 22 cm





Thermoplastic Dynamic Hand-Finger Splint

Ref : ORT-07

This is a hand and wrist orthosis made of thermoplastic material, covered with terry cloth inside. Keeps the wrist and the fingers in passive extension. Allows for active finger flexion. Provides early controlled motion after tendon repair in extensor tendon injury and prevent the formation of adhesions.



Size;

Left / Right Available.

M : 14 - 18 cm

L : 18 - 22 cm



Thermoplastic Kleinert Hand-Finger Splint

Ref : ORT-09

Made of thermoplastic material, covered with terry cloth inside. Allows passive flexion and active extension. Provides early controlled motion after tendon repair in flexor tendon injury and prevent the formation of adhesions.



Size;

Left / Right Available.

M : 14 - 18 cm

L : 18 - 22 cm



Finger Splints



Baseball Splint Frog Splint Finger Stabilizer Splint Straight Finger Splint

These are finger splints made of special aluminum, covered with plastizode inside. Provides finger immobilization in extension. Used for finger fractures, tendon and nerve injuries, strains and sprains.

Size;

Ref : ORT-02B

Ref : ORT-02K

Ref : ORT-02S

S : 4,5 - 5 cm

M : 5 - 6 cm

L : 6 - 7 cm

Ref : ORT-02D

Standard Size.

Ref: ORT-02B
Ref: ORT-02K
Ref: ORT-02S
Ref: ORT-02D



Mallet Finger Splint

Ref : ORT-01

This is a finger splint made of special plastic material, keeping the DIP joint in extension. The fingers are free after the PIP joint. Used for extensor tendon injuries, fractures, strains and sprains.

Size;

0 : 4 - 4,5 cm

1 : 4,5 - 5 cm

2 : 5 - 5,5 cm

3 : 5,5 - 6 cm

4 : 6 - 6,5 cm

5 : 6,5 - 7 cm

6 : 7 - 8 cm



Butonier Splint

Ref : ORT-03

This is a finger splint made of thermoplastic material and which keeps the finger in extension according to the 3-point principle.

Size;

S : 4,5 - 5 cm

M : 5 - 6 cm

L : 6 - 7 cm





Medical Corsets



Chest Corset with Steel Support

Ref: AO-23


This corset has the same specifications with AO-20 chest corset, the only difference being the steel underwire within the sternal area. Used for the same cases as Ref AO-20.

Indications;

- 1-Sternal support after coronary surgery,
- 2-Costa and sternum fractures,
- 3-After sternal surgery and trauma.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Chest Corset

Ref: AO-20


Manufactured with flexible, respirable and sweat-preventing cotton fabric. The sternal area is supported with terry cloth. It provides compression to the rib cage and supports the sternum.

Indications;

- 1-Sternal support after coronary surgery,
- 2-Costa and sternum fractures,
- 3-After sternal surgery and trauma.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Abdominal Corset

Ref: AO-25


This is the abdominal support corset manufactured with flexible, respirable and sweat-preventing cotton fabric.

Indications;

- 1-After abdominal surgery and birth,
- 2-Peritoneal injury,
- 3-Abdominal sagging.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm




Abdominal Corset (16 cm)

Ref: AO-31

This is the abdominal support corset manufactured with flexible, respirable and sweat-preventing cotton fabric. It differs from the AO-25 corset with its 16-cm width.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Woolen Corset

Ref: AO-26

Made of wool. Provides the waist area with warmth and compression

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Umbilical Hernia Corset


Ref: AO-24

Manufactured with flexible, respirable and sweat-preventing cotton fabric. It includes a pad designed especially for providing compression for the herniated area. Used for umbilical hernia.



Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Colostomy Corset (16 cm)


Ref: AO-29

Manufactured with flexible and sweat-preventing cotton fabric. Used in colostomy cases in order to protect the colostomy bag. It differs from AO-30 in that it is used after the removal of the stitches.



Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Colostomy Corset (26 cm)


Ref: AO-30

Manufactured with flexible and sweat-preventing cotton fabric. Used in colostomy cases in order to protect the colostomy bag. It is used after surgical operation.



Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Pregnancy Corset

Ref: AO-27

Manufactured with flexible and respirable cotton fabric, this is an abdominal support corset which might be used before and after giving birth. If used during pregnancy, it supports the abdomen and prevents lumbar pain and backaches. In daily life, it makes the patient more comfortable and confident while moving. After giving birth, corset is worn with the wider part in front; this prevents abdominal sagging and ensures that the patient returns to her daily activities easily and safely.



Size;

S : 85 - 95 cm
M : 95 - 105 cm
L : 105 - 115 cm
XL : 115 - 125 cm
XXL: 125 - 135 cm



Sacro-Iliac Corset

Ref: AO-75

Manufactured with flexible, respirable and sweat-preventing material. The perineal fasteners prevent the corset from shifting upwards.



Indications;

- 1-Injury or post-natal sacroiliac arthrosis or deviation of symphysis pubis,
- 2-Sacroiliac joint pain,
- 3-Post-operation after symphysis resection (for protection).

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Neoprene Lumbosacral Corset (with strap)

Ref: 201

Made of neoprene material. Supported with an extra outer strap. Provides warmth and compression.




Indications;

- 1-Lumbago,
- 2-Acute or chronic lumbar pain.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Lumbostad Corset (26 cm) with Strap

Ref: AO-60

Manufactured with flexible, respirable and sweat-preventing cotton fabric. Involves steel and flexible underwire. Its length is 26 cm. Relieves the pressure off the discs with the aid of transverse straps.




Indications;

- 1-Lumbago, degenerative disc lesions,
- 2-Spondylodesis, spondylolisthesis,
- 3-Acute or chronic lumbar pain.
- 4-Post-operation (for immobilization).

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Lumbosacral Corset (32 cm) with Strap

Ref: AO-80

Manufactured with respirable and sweat-preventing cotton fabric. Involves steel and flexible underwire. Relieves the pressure off the discs with the aid of transverse straps. It is 32 cm long and supports up to the waist and lower back.




Indications;

- 1-Lumbago, degenerative disc lesions,
- 2-Spondylodesis, spondylolisthesis,
- 3-Acute or chronic lumbar pain.
- 4-Post-operation (for immobilization).

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Lumbostad Corset With Pad

Ref: AO-63

Manufactured with respirable and flexible cotton fabric. Includes 4 pieces of steel underwire and pads made of sponge fabric in the rear lumbar area. Its length is 26 cm.




Indications;

- 1-Lumbago, degenerative disc lesions,
- 2-Spondylodesis, spondylolisthesis,
- 3-Acute or chronic lumbar pain.
- 4-Post-operation (for immobilization).

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Lumbosacral Corset with Special Support

Ref: AO-70

Manufactured with respirable and sweat-preventing cotton fabric. This is a two-piece corset with 4 pieces of inner underwire. On the exterior, it has an anatomical extra support made of lamination. It supports the lumbar area and, relieving the pressure on the discs, reduces pain.




Indications;

- 1-Lumbago, disc lesions,
- 2-Spondylodesis, spondylolisthesis,
- 3-Post-operation and -trauma (for protection and immobilization).

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Lumbostad Corset (26 cm) with Abdominal Support


Ref: AO-61

Manufactured with flexible, sweat-preventing cotton fabric. In the back, it has 4 pieces of anatomically convenient steel underwire. The front side is composed of two pieces. The lower and thicker piece supports the lower abdomen. The thinner piece on the front side is fastened in a transverse manner.



Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Woven Lumbosacral Corset

Ref: AO-85

It is manufactured with respirable woven fabric. In the back, it has 4 pieces of anatomically convenient steel underwire. It supports the lumbar area and relieves the pain.




Indications;

1-Lumbago, disc lesions,
2-Spondylodesis,
spondylolisthesis,
3-Post-operation and
-trauma (for protection
and immobilization).

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Woven Lumbosacral Corset with Extra Support

Ref: AO-86

It is manufactured with respirable woven fabric. In the back, it has 4 pieces of anatomically convenient steel underwire and plastic support. It supports the lumbar area and, relieving the pressure on the discs, reduces pain.




Indications;

- 1-Lumbago, disc lesions,
- 2-Spondylodesis, spondylolisthesis,
- 3-Post-operation and -trauma (for protection and immobilization).

Size;


S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Woven Lumbosacral Corset With Pad

Ref: AO-87

This is a corset made of flexible woven fabric, which is respirable and has a specially designed massage pad in the rear lumbar side. The massaging accelerated the blood flow. It stabilizes the sacroiliac and lumbosacral joints and relieves the painful area with the compression it applies, reducing the muscle tone.




Indications;

- 1-Symphysiolysis
- 2-Spondylosis, degenerative joint alterations,
- 3-Sacroiliac joint irritation,
- 4-Lumbar pain caused by intervertebral joint disorders.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm






Woven Lumbosacral Corset with Support


Ref: AO-88

Manufactured with breathing and flexible woven fabric. Includes 4 anatomically designed, inclined steel supports in the rear lumbar area.




Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm




Woven Lumbosacral Corset with Support, with Straps Ref: AO-89

Manufactured with breathing and flexible woven fabric. Includes 4 anatomically designed, inclined steel supports in the rear lumbar area. Differently, it has an extra strap which increases the stabilization of the lumbar area.



Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm






Dorsolumbar Corset (32 cm)

Ref: AO-90

Manufactured with breathing and sweat-preventing cotton fabric. It has steel underwire on both sides of the vertebral column extending up to the back. It has an extra strap with transverse fasteners in the lumbar area.




Indications;

- 1-Vertebral trauma,
- 2-Osteoporosis,
- 3-Mild kyphos.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Dorsolumbar Corset (26 cm)

Ref: AO-92

Manufactured with breathing and sweat-preventing cotton fabric. It has steel underwire on both sides of the vertebral column extending up to the back. It has an extra strap with transverse fasteners in the lumbar area. Since the lumbar part is 26 cm, short patients might also use it comfortably.




Indications;

- 1-Vertebral trauma,
- 2-Osteoporosis,
- 3-Mild kyphos.

Size;

S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





Dorsolumbar Corset (26 cm)

Ref: AO-91

Manufactured with breathing and sweat-preventing cotton fabric. It has steel underwire on both sides of the vertebral column extending up to the back. The lumbar part is 26 cm. Different from the common dorsolumbar corsets, the axillary strap is made of terry pads, which provide a more comfortable usage, and it consists of a buckle system for easier wear. In order to prevent axillary pressure, height adjustment is possible through the velcro fasteners in the dorsal area. The rear part of the corset is softened with terry pads. The cover is manufactured so as to allow for easy wear.




Indications;

- 1-Vertebral trauma,
- 2-Osteoporosis,
- 3-Mild kyphos.

Size;

- S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm



Posturex Bandage

Ref: AO-50

Manufactured with flexible, breathing and sweat-preventing cotton fabric. It has plastic underwire on both sides of the vertebral column. By preventing posture disorders, it provides sensory feedback.




Indications;

- 1-Posture disorders,
- 2-Weakness of back and shoulder muscles.

Size;

- S : 75 - 85 cm
M : 85 - 95 cm
L : 95 - 105 cm
XL : 105 - 115 cm
XXL: 115 - 125 cm





TLSO Corset

Ref: AO-95

This is a body orthosis made of sponge fabric with terry cloth inside. There is thermoplastic material, which provides a rigid support, in the abdominal and dorsal parts of the corset. Like the dorsolumbar corset, this one is also easy to wear, but what is different is that the support in the sternum restrict the flexion of the body movements. The support in the sternum can be lengthened and shortened in accordance with the patient's body.




Indications;

- 1-Vertebral fractures,
- 2-Post-operation,
- 3-Compression fractures.

Size;

- S : 60 - 85 cm
M : 85 - 110 cm
L : 110 - 140 cm



LSO Corset

Ref: AO-96

This is a lumbosacral orthosis made of sponge fabric with terry cloth inside. There is thermoplastic material, which provides a rigid support, in the abdominal and lumbar parts of the corset. It restricts the motion of the lumbar vertebrae to a certain extent.




Indications;

- 1-Post-operation in lumbar vertebra fractures,
- 2-Stable compression fractures
- 3-Post-trauma.

Size;

- S : 60 - 85 cm
M : 85 - 110 cm
L : 110 - 140 cm





Hyperextension Corset (4-point)

Ref: ORT-12

This is a body orthosis made of special aluminum and its parts which are in contact with the skin are covered with laminated pads. Due to the screw system installed on the orthosis, length and width adjustment is possible. It stabilizes the hyperextension of the vertebra. Its rotation is restricted with subclavicular extension pads. It restricts lateral flexion with stable skeletal structure. Comes in 2 sizes.

Indications;

- 1-Stable compression fracture of lumbar and lower thoracic vertebrae,
- 2-Immobilization after intervertebral disc surgery,
- 3-Axial stabilization after osteoporotic compression fracture of thoracic and lumbar vertebrae,
- 4-After non-stable vertebral body fracture surgery,
- 5-Non-operable spinal tumor and its metastases,
- 6-After the surgery for vertebral body tumor and its metastases, for post-operation orthotic care and orthosis,
- 7-After spinal decompression surgery where internal fixation was performed or not.

Size;

- S : 40 - 49 cm
L : 50 - 60 cm





Hyperextension Corset (3-point)

Ref: ORT-11

This is a body orthosis made of special aluminum and its parts which are in contact with the skin are covered with laminated pads. Due to the screw system installed on the orthosis, length and width adjustment is possible. It restricts the vertebral movement to a certain extent. Comes in 2 sizes.



Indications;

- 1-Stable compression fracture of lumbar, middle and lower thoracic vertebrae,
- 2-Hyperkyphos in Scheuermann's disease,
- 3-For the conservative treatment of non-operable vertebral body fractures in the lumbar and thoracic vertebrae,
- 4-Temporary post-operation treatment.

Size;

S : 40 - 49 cm
L : 50 - 60 cm



Hyperextension Corset

Ref: ORT-13

It is made of a special kind of aluminum; the parts which are in contact with the skin are covered with terry pads; in accordance with the 3-point principle, this orthosis restricts the flexion of the body and encourages the hyperextension position. It restricts the vertebral movements to a certain extent. Due to the screw system installed on the orthosis, length and width adjustment is possible. Comes in 2 sizes.



Indications;

- 1-Vertebral fractures and trauma,
- 2-Compression fractures in osteoporosis,
- 3-Kyphos,
- 4-Post-operation..

Size;

S : 40 - 49 cm
L : 50 - 60 cm





Inguinal Truss (Single Sided)

Ref: K-100

Manufactured with flexible cotton fabric. There is a special pad applying pressure on the inguinal area. This is used in cases of single sided inguinal hernia. The buckle system makes it easy to use.




Size;

Left / Right Available.

S : 75 - 90 cm

M : 90 - 105 cm

L : 105 - 120 cm



Inguinal Truss (Double Sided)

Ref: K-120

Manufactured with flexible cotton fabric. There is a special pad applying pressure on the inguinal area. This is used in cases of double sided inguinal hernia. The buckle system makes it easy to use.



Size;

S : 75 - 90 cm

M : 90 - 105 cm

L : 105 - 120 cm



Child Umbilical Hernia Truss

Ref: ÇG-121

Manufactured with flexible fabric. It has a pad applying compression to the umbilical hernia. It is of standard size but adjustable. Used for umbilical hernia in children.



Size;


Standard Size.



Child Inguinal Hernia Truss (Single Sided)

Ref: ÇK-101

Manufactured with flexible fabric. It has a pad applying compression to the hernia. Perineal strap prevents sliding. It is of standard size but adjustable. It has a left and a right side. Used for inguinal hernia in children.



Size;

Standard Size.

Left / Right Available.





Lower extremity



Patella Knee Support

Ref: 101

Manufactured anatomically with neoprene material. Perimeter of the patella is open and pad-supported. Provides ease of use with the Velcro fasteners in the upper calf.

Indications;

- 1-Arthritis,
- 2-Increase of intra articular liquid,
- 3-Pain induced by knee cap,
- 4-Post-operation and -trauma (for protection and rehabilitation),
- 5-Protection during sports activities,
- 6-For heat and compression.

Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL : 44 - 48 cm



Patella and Ligament Knee Support

Ref: 102

Manufactured anatomically with neoprene material. Perimeter of the patella is open and pad-supported. The flexible underwire on the sides increase stabilization and thus decrease the load on the knee. Provides ease of use with the Velcro fasteners in the upper calf.

Indications;

- 1-Arthritis, bursitis,
- 2-Increase of intra articular liquid,
- 3-Pain induced by knee cap,
- 4-Meniscus problems,
- 5-1st degree twist and ligament instability,
- 6-Post-operation (for protection).

Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL : 44 - 48 cm



Simple Knee Support

Ref: 100

Manufactured anatomically with neoprene material. Provides ease of use with the Velcro fasteners in the upper calf. Supports the knee joint. Reduces pain by decreasing the load on the knee.

Indications;

- 1-Arthritis,
- 2-Increase of intra articular liquid,
- 3-Post-operation and -trauma (for protection and rehabilitation),
- 4-Protection during sports activities,
- 5-For heat and compression.

Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL : 44 - 48 cm





Patella and Ligament Knee Support with Velcro

Ref: 103

Manufactured anatomically with neoprene material. This has the same specifications as Ref 102. Differently, there is more stabilization due to its longer size and by the help of Velcro fasteners in the upper and lower calf.

Indications;

- 1-Arthritis, bursitis, tendinitis,
- 2-Sprains and Strains
- 3-Meniscus problems,
- 4-Post-operation (for protection).
- 5-Ligament instability,
- 6-Sports injuries.

Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL : 44 - 48 cm



Hinged Knee Support (Open Top)

Ref: 104

Manufactured anatomically with neoprene material. Aluminum joints on both sides provide maximum support and allow for controlled motion. Prevents knee hyperextension. Provides ease of use with the opening from the front side in the upper calf. There are Velcro fasteners on the upper and lower calves increasing stabilization.

Indications;

- 1-Ligament injury and instability,
- 2-Meniscus lesions,
- 3-1st or 2nd degree sprains and strains
- 4-Post-operation and -trauma (for protection and immobilization).

Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL : 44 - 48 cm



Hinged Knee Support (Open Front)

Ref: 104-A

This has the same specifications as Ref 104. Since can be completely opened from the front side, it provides ease of use in overweight and elderly patients. The presence of a rear popliteal window folds while bending the knee, preventing the disturbance and pressure. Used for the same cases as Ref104.



Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL : 44 - 48 cm






Patellar Stabilization Knee Support

Ref: 107

Made of neoprene material. There is flexible underwire on both sides. Moreover, the silicone pad around the patella support the patella from the inside or outside according to the situation. It has a left and a right side.



Indications;

- 1-Patella instability,
- 2-Patellar medialization,
- 3-Patellar lateralization,
- 4-Knee cap cartilage softening.

Size;

Left / Right Available.

S : 31 - 34 cm

M : 34 - 37 cm

L : 37 - 40 cm

XL : 40 - 44 cm

XXL: 44 - 48 cm



Knee Support With Pad

Ref: 105

Manufactured anatomically with neoprene material. Patella and its perimeter is supported with an especially designed laminated pad.

Indications;

- 1-Sportif aktiviteler sırasında koruma amaçlı,
- 2-Patellayı ve çevre yumuşak dokuları travmaya karşı koruma amaçlı,
- 3-Dizdeki bursitleri koruma amaçlı kullanılır.

Size;

S : 31 - 34 cm

M : 34 - 37 cm

L : 37 - 40 cm

XL : 40 - 44 cm

XXL: 44 - 48 cm



Knee Support for Post-arthroscopy

Ref: 106

Manufactured anatomically with neoprene material. The inner pocket system is designed to place hot and cold application gel. Since can be completely opened from the front side, it provides ease of use in post-operation and elderly patients. Threesome Velcro fasteners increase stabilization.

Indications;

- 1-Treatment of post-arthroscopy edema and pain,
- 2-Any condition where hot or cold application is preferred.

Size;

S : 31 - 34 cm

M : 34 - 37 cm

L : 37 - 40 cm

XL : 40 - 44 cm

XXL: 44 - 48 cm





Knee Support for Ant. Cruciate Ligament

Ref : 109

Manufactured anatomically with neoprene material. Perimeter of the patella is open and pad-supported. There is flexible underwire on both sides. There are transverse straps adhering with 3-point principle. Provides stabilization for front transverse ligaments and patella

Indications;

- 1-Transverse ligament injury and knee instability,
- 2-Patellar stabilization,
- 3-Meniscus lesions,
- 4-Post-operation and -trauma (for protection and stabilization).

Size;

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm



Patellar Tendon Band

Ref : 110


Manufactured anatomically with neoprene material. There is a special pad supporting the patellar tendon. Does not restrict joint motion. The compression effect reduces tension and pain in the tendon.

Indications;

- 1-Osgood schlatter's disease,
- 2-Patellar tendon disorders.

Size;

Standard Size.



Knee Support for Patellar Tendon

Ref : 111


Manufactured anatomically with neoprene material. There is a special pad supporting the patellar tendon. The compression applied reduces pain.

Indications;

- 1-Osgood schlatter's disease,
- 2-Patellar tendon disorders.

Size;

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm





Woolen Knee Support

Ref: 112

Made of wool. Provides warmth and compression

Size;

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm



Woven Knee Support

Ref: 113

Manufactured with flexible and sweat-preventing woven fabric.

Indications;

- 1-Arthritis,
- 2-Increase of intra articular liquid,
- 3-Post-operation (for protection).
- 4-Protection during sports activities.

Size;

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm






Woven Patella Knee Support

Ref: 119

Manufactured with flexible and sweat-preventing woven fabric, anatomically convenient for the shape of the knee. The perimeter of the patella is supported with a silicone pad.



Indications;

- 1-Osteoarthritis,
- 2-Increase of intra articular liquid,
- 3-Pain caused by patella,
- 4-Post-operation and -trauma (for protection and rehabilitation).

Size;


S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm



Patella and Single Ligament Woven Knee Support

Ref : 120

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. The perimeter of the patella is supported with silicone. There is a single underwire on both sides. Supports the knee joint, reduces pain by relieving the load on the joint.



Size;

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm






Patella and Ligament Woven Knee Support

Ref : 114

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. The perimeter of the patella is supported with silicone. There is flexible underwire on both sides. Supports the knee, reduces pain by relieving the load on the joint.



Indications;

- 1-Osteoarthritis, bursitis,
- 2-Meniscus lesions,
- 3-Increase of intra articular liquid,
- 4-Pain induced by knee cap,
- 5-Ligament instability,
- 6-Post-operation and -trauma (for protection and stabilization).

Size;


- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm



Woven Ligament Knee Support

Ref: 118

Manufactured with flexible and sweat-preventing woven fabric, anatomically convenient for the shape of the knee. There are two ligament supports on both sides made of flexible spiral underwire



Indications;

- 1-Osteoarthritis,
- 2-Muscle and ligament rupture,
- 2-Bursitis,
- 4-Post-operation (for protection).

Size;

- S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm






Patellar Stabilization Woven Knee Support

Ref: 115

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. There is flexible underwire on both sides. The silicone pad around the patella support the patella from the inside or outside according to the situation



Indications;

- 1-Patella instability,
- 2-Patellar medialization,
- 3-Patellar lateralization,
- 4-Chondromalacia patella cases.

Size;

Left / Right Available.

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm



Woven Hinged Knee Support

Ref: 116

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. The perimeter of the patella is supported with silicone and there are aluminum joints on both sides. There are Velcro fasteners on the upper and lower calves increasing stabilization.



Indications;

- 1-Ligament injury,
- 2-Meniscus rupture,
- 3-Arthritis, bursitis, tendinitis,
- 4-Post-operation and -trauma (for protection and support).

Size;

S : 31 - 34 cm
M : 34 - 37 cm
L : 37 - 40 cm
XL : 40 - 44 cm
XXL: 44 - 48 cm





Aurastrap (Hinged)

Ref: 170

This is a functional hinged rigid brace allowing for flexion and extension limitation, designed according to the shape of the knee and manufactured with light aluminum. It provides the most convenient stabilization with the 4-point principle. It prevents tibial displacement, also preventing knee hyperextension.




Indications;

- 1-Isolated rupture of ant cruciate ligament or combined ant and post cruciate ligament ruptures ,
- 2-Instability caused by collateral ligament rupture,
- 3-Conservative treatment of combined instabilities.

Size;

Left / Right Available.


Perimeter of the Part below 16 cm from the Knee	Perimeter of the Part above 16 cm from the Knee
S : 34 - 38 cm	42 - 48 cm
M : 38 - 42 cm	48 - 54 cm
L : 42 - 46 cm	54 - 60 cm
XL : 46 - 50 cm	60 - 67 cm
XXL: 50 - 54 cm	67 - 74 cm



Aurastrap (Angle Adjusted)

Ref: 171

This is a functional angle adjusted rigid brace allowing for flexion and extension limitation, designed according to the shape of the knee and manufactured with light aluminum. Flexion (0-15-30-60-90-120) and extension (0-5-15-30) angle adjustment is possible. It provides the most convenient stabilization with the 4-point principle. It prevents tibial displacement, also preventing knee hyperextension.




Indications;

- 1-Isolated rupture of ant cruciate ligament or combined ant and post cruciate ligament ruptures ,
- 2-Instability caused by collateral ligament rupture,
- 3-Conservative treatment of combined instabilities.

Size;

Left / Right Available.

Perimeter of the Part below 16 cm from the Knee	Perimeter of the Part above 16 cm from the Knee
S : 34 - 38 cm	42 - 48 cm
M : 38 - 42 cm	48 - 54 cm
L : 42 - 46 cm	54 - 60 cm
XL : 46 - 50 cm	60 - 67 cm
XXL: 50 - 54 cm	67 - 74 cm





Immobilizer

Ref: AO-45/ AO-55/ AO-65

Manufactured with a sponge fabric with terry cloth inside. Provides knee immobilization in extension. Designed anatomically, it has aluminum supports providing stabilization in the posterior, medial and lateral. Reinforced with Velcro fasteners. It comes in standard size, but available in 45 cm, 55 cm, 65 cm lengths.




Indications;

- 1-Helpful for post-operation recovery, protection and immobilization.
- 2-For recovery after ligament injury,
- 3-Post-operation (for meniscus injury).

Size;

AO-45 : 45 cm
AO-55 : 55 cm
AO-65 : 65 cm



Immobilizer (20 degree)

Ref: AO-47/ AO-57/ AO-67

Manufactured with a sponge fabric with terry cloth inside. Provides the knee with immobilization at a 20 degree flexion. Designed anatomically, it has aluminum supports providing stabilization in the posterior, medial and lateral. Reinforced with Velcro fasteners. It comes in standard size, but available in 45 cm, 55 cm, 65 cm lengths.



Size;

AO-47 : 45 cm
AO-57 : 55 cm
AO-67 : 65 cm





Angle-Adjusted Brace

Ref: 740

This is a degree-adjustable knee orthosis made of an aluminum bar and fabric with terry cloth interior which contacts the skin. Flexion (0-15-30-60-90-120) and extension (0-5-15-30) angle adjustment is possible. Comes in 2 sizes.



Indications;

- 1-Gradual mobilization after operation,
- 2-Controlled and gradual joint mobilization after ligament injury operations,
- 3-Help for ambulation.

Size;

SHORT : 47 cm
LONG : 58 cm



Angle-Adjusted Brace

Ref: 745

This is a degree-adjustable knee orthosis made of an aluminum bar and fabric with terry cloth interior which contacts the skin. Flexion (0-15-30-45-60-75-90-105-120) and extension (-0-15-15-30-45-60-75-90) angle adjustment is possible. Comes in 2 sizes.



Indications;

- 1-Gradual mobilization after operation,
- 2-Controlled and gradual joint mobilization after ligament injury operations,
- 3-Help for ambulation.

Size;

SHORT : 47 cm
LONG : 58 cm





Thigh Support

Ref: 500

Manufactured anatomically with neoprene material. Provides ease of use with the Velcro fasteners in the upper part. Supports the upper calf. It accelerates the blood flow and reduces pain with heat and compression.

Indications;

- 1-Muscle injuries to the thigh,
- 2-Thigh pain,
- 3-Protection during and after sports activities,
- 4-For providing heat and compression.

Size;

- S : 44 - 48 cm
M : 48 - 52 cm
L : 52 - 56 cm
XL : 56 - 62 cm
XXL: 62 - 68 cm



Calf Support

Ref: 501

Manufactured anatomically with neoprene material. Provides ease of use with the Velcro fasteners in the upper part. Supports the lower calf. It accelerates the blood flow and reduces pain with heat and compression.

Indications;

- 1-Muscle injuries to the thigh,
- 2-Calf pain,
- 3-Protection during and after sports activities,
- 4-For providing heat and compression.

Size;

- S : 29 - 32 cm
M : 32 - 35 cm
L : 35 - 39 cm
XL : 39 - 43 cm
XXL: 43 - 48 cm





Hip Support

Ref: 502

Manufactured anatomically with neoprene material. Provides compression around the thigh and hip joint circumference and protects against trauma. Ease of use with the Velcro fastener on the end section of thigh.

Indications;

- 1-Muscle injuries in thigh and hip circumference,
- 2-Protection of the soft tissue in the hip circumference,
- 3-Hip circumference bursitis,
- 4-For providing heat and compression.

Size;

Left / Right Available.

S : 44 - 48 cm

M : 48 - 52 cm

L : 52 - 56 cm

XL : 56 - 62 cm

XXL: 62 - 68 cm



Slimming Tights

Ref: 503

Made of neoprene material. Ease of use with the zipper.

Indications;

- 1-For providing heat and compression,
- 2-For inguinal problems,
- 3-For slimming around the waist and hip line.

Size;

S : 75 - 85 cm

M : 85 - 95 cm

L : 95 -105 cm

XL :105 -115 cm

XXL:115 -125 cm





Ligament Ankle Support

Ref: 404

Made of neoprene material. Around the malleolus, the stabilizing inner and outer plastic plates are produced separately in accordance with the position of the malleolus. It can be completely opened, thus it is easy to use. The figure-8 elastic transverse strap provides extra support. It supports the ankle, reduces edema and pain by compression.



Indications;

- 1-Ligament instability and rupture,
- 2-Protection and edema prevention,
- 3-Post-operation and –injury
(for protection and support).

Size;

S : 21 - 23 cm
M : 23 - 25 cm
L : 25 - 27 cm
XL : 27 - 29 cm
XXL: 29 - 31 cm



Simple Ankle Support

Ref: 401

Manufactured anatomically with neoprene material. Ease of use with the Velcro fasteners around the ankle. Supports the ankle, thus reduces pain.

Indications;

- 1-1st degree twist and tonus,
- 2-Small fractures,
- 3-Compression in edema after injury,
- 4-Protection and stabilization after twisting.

Size;

S : 21 - 23 cm
M : 23 - 25 cm
L : 25 - 27 cm
XL : 27 - 29 cm
XXL: 29 - 31 cm



Cross Strap Ankle Support

Ref: 402

Manufactured anatomically with neoprene material. The figure-8 elastic transverse strap provides extra stabilization and compression. Used for the same cases as Ref401.

Size;

S : 21 - 23 cm
M : 23 - 25 cm
L : 25 - 27 cm
XL : 27 - 29 cm
XXL: 29 - 31 cm





Dorsi-Flexion Bandage

Ref: 411

Made of neoprene material. This is an ankle band used for fallen arch deformity, which keeps the ankle within the desired stretch in the dorsal flexion.



Size;


Standard Size.



Malleol Ankle Support

Ref: 405

Made of neoprene material. The inner and outer malleolus is supported with anatomically designed pads. Supported with the elastic band wrapped in the shape of an eight Supports the ankle, relieves the edema and thus reduces pain.



Indications;

- 1-Mild ligament injury and instability,
- 2-Osteoarthritis and arthritis,
- 3-Post-operation and –injury
(for protection and support).

Size;

- S : 21 - 23 cm
- M : 23 - 25 cm
- L : 25 - 27 cm
- XL : 27 - 29 cm
- XXL: 29 - 31 cm



Eight Bandage Ankle Support

Ref: 406

Made of respirable, soft cotton fabric, it can comfortably be worn inside the shoe. Supports the ankle, relieves the edema.

Indications;

- 1-Helps ligament instability,
- 2-Ankle pain.

Size;

- S : 21 - 23 cm
- M : 23 - 25 cm
- L : 25 - 27 cm
- XL : 27 - 29 cm
- XXL: 29 - 31 cm





Woven Ankle Support

Ref: 407

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. The compression effect reduces edema and pain.

Indications;

- 1-1st degree strain and sprain.
- 2-After minor fractures,
- 3-Edema after injury,
- 4-For protection after sprains,
- 5-Post-operation (for protection).

Size;

S : 19 - 21 cm
M : 21 - 23 cm
L : 23 - 25 cm
XL : 25 - 27 cm
XXL: 27 - 29 cm



Woven Varus Ankle Support

Ref: 403

Manufactured with flexible and sweat-preventing woven fabric. With the elastic bandage, it treats the varus position (stepping with the ankle bent outwards) passively.

Size;

S : 19 - 21 cm
M : 21 - 23 cm
L : 23 - 25 cm
XL : 25 - 27 cm
XXL: 27 - 29 cm






Achilles Tendon Ankle Support

Ref: 409

Designed anatomically, manufactured with flexible and sweat-preventing woven fabric. There are silicone supports on both sides of the Achilles' tendon and under the heel. Also, there is a portable heel cup placed under the other heel in order to create evenness between two legs. Protects the heel and the Achilles' tendon, reduces the tension in the tendon.



Indications;

- 1-Aşıl tendon akut ve kronik iltihaplanmalarında,
- 2-Aşıl tendon ağrısı,
- 3-Retrokalkaneal bursit,
- 4-Ameliyat veya yaralanma sonrası kullanılır.

Size;


- S : 19 - 21 cm
M : 21 - 23 cm
L : 23 - 25 cm
XL : 25 - 27 cm
XXL: 27 - 29 cm



Woven Malleol Ankle Support

Ref: 408

Manufactured with flexible and sweat-preventing woven fabric, it has anatomically designed silicone pads around malleolus. Supports the ankle, reduces pain and edema with compression.



Indications;

- 1-Mild ligament injury and instability,
- 2-Osteoarthritis and arthritis,
- 3-Post-operation and –injury (for protection and support).

Size;

- Has Left-Right Sides.
S : 19 - 21 cm
M : 21 - 23 cm
L : 23 - 25 cm
XL : 25 - 27 cm
XXL: 27 - 29 cm






Lace-up Ankle Support

Ref: 417

Manufactured with respirable perforated material. There is a ligament support on both sides. The lace system allows a complete wrapping of the ankle, thus stabilization. Supported with an extra figure-8 bandage. The rear part is zipped for easier use.



Indications;

- 1-Pre- and post-operation for ligament injury,
- 2-For compression and support after sprain and injury,
- 3-For protection in sports injuries.

Size;

- S : 21 - 23 cm
M : 23 - 25 cm
L : 25 - 27 cm
XL : 27 - 29 cm
XXL: 29 - 31 cm



Stabilization Ankle Brace

Ref: 412

The exterior is made of special rigid plastic, while the interior is neoprene. Provides the ankle with inner and outer stabilization. Comes in standard size, adjustable with Velcro fasteners around the ankle.

Indications;

- 1-Chronic ankle instability,
- 2-Pre- and post-operation for ligament injury,
- 3-For compression and support after twisting,
- 4-For protection in sports injuries.

Size;

Standard Size.





Ankle Brace with Gel

Ref: 410

Made of plastic material anatomically designed on the outer surface, this is an ankle orthosis whose interior enables cold treatment with gel pads. Provides the ankle with medial and lateral stabilization.



Indications;

- 1-Chronic ankle instability,
- 2-Pre- and post-operation for ligament injury,
- 3-For compression and support after sprain,
- 4-For protection in sports injuries,
- 5-Stabilization after taking out the plaster of a fracture.

Size;

Standard Size.



Ankle Brace With Pad

Ref: 416

Made of plastic material anatomically designed on the outer surface, this is an ankle orthosis whose interior is made of sponge pads. Provides the ankle with inner and outer stabilization.



Size;

Standard Size.





R.O.M. Walker

Ref: 413

Manufactured with a sponge fabric with terry cloth inside. This is an ankle orthosis which helps early mobilization, dorsal and plantar flexion degree is adjustable. By the roomwalker, it is provided a plantar flexion of ankle in degrees 0-10-20-30-40 and dorsiflexion of ankle in degree 20-10-0 in each angle.



Indications;

- 1-Tendon, ligament and soft tissue injury (for post-operation stabilization)
- 2-Chronic ligament instability,
- 3-Stable distal tibia and fibula fractures or distal tibia and fibula fractures where an internal fixer is used,
- 4-Achilles' tendon rupture,
- 5-Stable foot and ankle fracture.

Size;

- S : 34 - 37 number
M : 38 - 41 number
L : 42 - 45 number



Achilles Boot

Ref: 414

Manufactured with a sponge fabric with terry cloth inside. This is an ankle orthosis helpful for early mobilization, plantar flexion height adjustable by the aid of the pads which is three degree in the heel.

Indications;

- 1-Tendon, ligament and soft tissue injury (for post-operation stabilization)
- 2-Chronic ligament instability,
- 3-Stable distal tibia and fibula fractures or distal tibia and fibula fractures where an internal fixer is used,
- 4-Achilles' tendon rupture,
- 5-Stable foot and ankle fracture.

Size;

- S : 34 - 37 number
M : 38 - 41 number
L : 42 - 45 number





Foot Guard

Ref: ORT-15

This is a T strapped ankle orthosis made of thermoplastic material, the interior of which is covered with terry cloth. Ensuring that the ankle is fixed at 90 degrees, prevents rotation of the ankle with the T strap. Preferable for inpatients.

Size;

M : 35 - 39 number

L : 40 - 44 number



Refleks Afo

Ref: ORT-16

This is an ankle orthosis made of thermoplastic material. Positions the ankle at 90 degrees. Ease of use with the Velcro fastener going around the ankle. Suitable for use within shoes.

Indications;

- 1-Fallen arch deformity caused by nerve injuries such as hemiplegia, lumbar disc hernia.
- 2-Used in every situation where desired ankle angle is 90 degrees.

Size;

Left / Right Available.

S : 34 - 37 number

M : 38 - 41 number

L : 42 - 45 number





Hip Brace

Ref: 750

This is a hip abduction orthosis providing hip stabilization, anatomically designed and made of thermoplastic material. The interior which contacts the skin is covered with terry cloth. Flexion adjustment with 15-degree grades until 120 degrees is possible. Abduction angle is adjustable.




Indications;

- 1-Stabilization after total hip prosthesis or revision,
- 2-Femur head resection,
- 3-After femoral or acetabular osteotomy,
- 4-Cases of recurrent hip dislocation.

Size;

Has Left-Right Sides.	
Upper calf diameter	Waist diameter
S : 40 - 45 cm	60 - 75 cm
M : 45 - 52 cm	75 - 90 cm
L : 52 - 60 cm	90 - 105 cm
XL : 60 - 70 cm	105 - 120 cm





Pavlic Bandage

Ref: 760

Manufactured with polyamide material, the bootie, shoulder and chest circumference is softened with terry cloth. Positions the hip joint at a desired tension in flexion and abduction. Used for hip dislocation in children.

Size;

S : 33 - 41 cm
M : 41 - 51 cm
L : 51 - 61 cm



Frejka Orthosis

Ref: 770

Manufactured with 12-mm plastizode. Positions the hip joint in flexion and abduction. Used for hip dislocation in children.

Size;

S (1) : 27 - 32 cm
M (2) : 32 - 37 cm
L (3) : 37 - 42 cm
XL (4) : 42 - 47 cm
XXL (5) : 47 - 52 cm



Hip Abduction Orthosis

aurafix®
ORTHOPEDIC PRODUCTS



Denis Browne Splint

Ref: 780

This is a hip abduction orthosis made of aluminum. A pair of boots which are manufactured especially for the orthosis are fitted between both feet. The channel system is length-adjustable. Used for hip dislocation and pes equinovarus deformity.



Size;


Kısa : 29 - 41 cm
Uzun : 42 - 68 cm



Denis Browne Splint (Mobile)

Ref: 782

This is a more dynamic orthosis than Ref 780. The child can perform independent movement of the legs such as kicking or crawling. Prevents the risk of the boot sliding off of the foot and regional rash. The child is more active and sleeps more comfortably.



Size;

Kısa : 29 - 41 cm
Uzun : 42 - 68 cm



Denis Browne Boot

Ref: 785

This is the boot which can be fitted onto the Denis Browne splint; the bottom is stout leather and the top is leather with laces, toe and heel parts are open.

Size;

From size 12 to 23.





Ref: 3101
Standard Patella
Knee Support



Ref: 3103
Standard Patella and Ligament
Knee Support



Ref: 3104
Standard Hinged Knee Support



Ref: 3202
Standard Lumbostad Corset



Ref: 3607
Standard Wrist Splint



Ref: 3400
Standard Ligament
Ankle Support



Ref: 3608
Standard Wrist Splint with
Thumb Support



Ref: 3401
Standard Ankle Support



Elastic Ankle Support

Ref: AV-140

It s made of an elastic material . Elastic supports protect the joints and reduces injury risks.It reduces the pain with compression effect by increasing blood circulation.

Indications;

- 1- First degree of twisting and stretching up
- 2- After inquiry compression for edema
- 3- Protection after injuries

Size;

XS : 18 - 20 cm
S : 20 - 23 cm
M : 23 - 26 cm
L : 26 - 29 cm
XL : 29 - 32 cm



Elastic Calf Support

Ref: AV-141

It s made of an elastic material. Supports lower calf. It reduces the pain with compression effect by increasing blood circulation.

Indications;

- 1-Muscle tensions
- 2-Inflamation in upper calf
- 3-Protection during sports activities

Size;

XS : 29 - 32 cm
S : 32 - 35 cm
M : 35 - 38 cm
L : 38 - 41 cm
XL : 41 - 45 cm



Elastic Knee Support

Ref: AV-142

It s made of an elastic material .Protects knee joint . Reduces the pain with compression effect by increasing blood circulation.

Indications;

- 1-Muscle tensions
- 2-After inquiry compression for edema
- 3-Protection during sports activities

Size;

XS : 30 - 34 cm
S : 34 - 38 cm
M : 38 - 42 cm
L : 42 - 47 cm
XL : 47 - 53 cm





Foot Health Products



Silicone Medial & lateral Wedge (Adhesive)

Ref: 803

It is made of silicone material and the upper surface is covered with fabric. In accordance with the problem, it is placed on the inner or outer side of the shoe. Due to its adhesive nature, it can be fixed within the shoe.

Indications;

- 1-Calcanus varus,
- 2-Calcanus valgus,
- 3-O-Leg, X-Leg,
- 4-Meniscus degeneration.

Size;

- S : 34 - 36 number
M : 37 - 39 number
L : 40 - 42 number
XL : 43 - 45 number



Metatarsal Arch Support (Adhesive)

Ref: 800

It is made of silicone material and the upper surface is covered with fabric. It is placed within the proximal of the metatarsus heads. It relieves the pressure of the metatarsus heads and reduces pain. Due to its adhesive nature, it can be fixed within the shoe.

Indications;

- 1-Over-sensitivity of metatarsus head,
- 2-Metatarsalgia,
- 3-Callus developing under the metatarsus,
- 4-Morton's foot- hallux rigidus-claw toe deformity,
- 5-Acceptical necrosis of metatarsus head (Freiberg's disease).

Size;

- S : 34 - 39 number
M : 40 - 45 number



Longitudinal Arch Support (Adhesive)

Ref: 801

It is made of silicone material and the upper surface is covered with fabric. It supports the longitudinal arch and reduces the increased tension in the sole. It starts from the calcaneus' inner tubercle and continues until the rear of the metatarsus heads. Due to its adhesive nature, it can be fixed within the shoe.

Indications;

- 1-Pes planus,
- 2-Pes cavus,
- 3-When the longitudinal arch needs to be supported.

Size;

- S : 34 - 39 number
M : 40 - 45 number



Heel Cup (Adhesive)

Ref: 802

Üst yüzeyi kumaş kaplı, silikon malzemeden üretilmiştir. Kalkaneal epin olan bölge yumuşatılmıştır. Bu şekilde ağırlı noktaya yük binmesi engellenmiş olur. Yapışkan özelliği sayesinde ayakkabı içinde sabitlenebilir.

Indications;

- 1-Heel spur,
- 2-Plantar fasciitis,
- 3-Avescular necrosis of calcaneus (sever's disease),
- 4-Heel pain.

Size;

- S : 34 - 36 number
M : 37 - 39 number
L : 40 - 42 number
XL : 43 - 45 number





Silicone Insole Covered with Fabric

Ref: 804

It is made of silicone material and the upper surface is covered with fabric. It distributes the load on the arch of the foot evenly and relieve the pressure on the joints. Being fabric covered, it is easy to use. Used for the same cases as Ref-805.

Size;

XS : 35 - 36 number
S : 37 - 38 number
M : 39 - 40 number
L : 41 - 42 number
XL : 43 - 44 number
XXL : 45 - 46 number



Silicone Insole

Ref: 805

Made of silicone material. It distributes the standing body weight evenly and reduces load and stress on the joints. It absorbs the reaction force of the floor.

Indications;

- 1-Pes planus,
- 2-Plantar fasciitis,
- 3-Heel pain,
- 4-Metatarsus pain.

Size;

XS : 35 - 36 number
S : 37 - 38 number
M : 39 - 40 number
L : 41 - 42 number
XL : 43 - 44 number
XXL : 45 - 46 number



Silicone Heel Cup

Ref: 806

Made of silicone material. The area with the calcaneal epin is softened. It lessens the load on the heels which causes pain.

Indications;

- 1-Heel spur,
- 2-Plantar fasciitis,
- 3-Avascular necrosis of calcaneus (sever's disease),
- 4-Heel pain.

Size;

S : 34 - 37 number
M : 38 - 41 number
L : 42 - 45 number





Silicone Heel Cup (Optional)

Ref: 823

Made of silicone material. It has portable silicone cushions with 3 different hardness level for the calcaneal epin area. It lessens the load on the heels which causes pain. There are 3 hardness levels: Red (hardness), green (medium), blue (soft).

Indications;

- 1-Heel spur,
- 2-Plantar fasciitis,
- 3-Avescular necrosis of calcaneus (sever's disease),
- 4-Heel pain.

Size;

- S : 34 - 37 number
M : 38 - 41 number
L : 42 - 45 number



Starflex Insole

Ref: 807

This is an arch support made of special thermoplastic material and which has a leather-covered longitudinal and metatarsal support.

Indications;

- 1-Pes planus,
- 2-Pes kavus,
- 3-Transvers arkin desteklenmesi gereken durumlarda kullanılır.

Size;

Between sizes 35 - 44



Starflex insole (with Epine)

Ref: 808

It has the same specifications as Ref 807, the only difference being the softened epin area. In addition to the usage of Ref 807, it is used for cases of heel pain and heel spur.

Size;

Between sizes 35 - 44





Steel Insole

Ref: 811

This insole, manufactured with stout leather and supported with stainless steel bridge, has longitudinal and metatarsal support.

Indications;

- 1-Pes planus,
- 2-Pes cavus,
- 3-When the longitudinal arch needs to be supported.

Size;

Between sizes 34 - 46



Cork Insole

Ref: 809

This is an insole made of special cork and rubber mixture which has longitudinal and metatarsal support.

Indications;

- 1-Pes planus,
- 2-Pes cavus,
- 3-When the transverse arch needs to be supported.

Size;

Between sizes 35 - 44



Cork Insole (with Epine)

Ref: 810

It has the same specifications as Ref 809, the only difference being the softened epine area. In addition to the usage of Ref 809, it is used for cases of heel pain and heel spur.

Size;

Between sizes 35 - 44





Elastic Hallux Valgus Splint

Ref: AO-09

Made of elastic material. Preventing deviation of the toe to the lateral, it can provide a pulling force of the desired strength. Uses for hallux valgus deformity.

Size;

Standard Size.



Hallux Valgus Night Splint

Ref: AO-10

Manufactured with a sponge fabric with terry cloth inside. It prevents deviation of the toe to the lateral. Convenient for night use in hallux valgus deformity.

Size;

Standard Size.

Left / Right Available.



Silicone Hammer Toe Cushion

Ref: 813

Made of silicone material. The hammer is worn on the deformed toe like a cover. It applies pressure on the PIP joint and affects the passive treatment of the deformity.

Size;

Standard Size.

Left / Right Available.



Silicone Toe Spreader

Ref: 812

Made of silicone material. Suitable for day use. It is used for preventing the deviation of the toe to the lateral in hallux valgus deformity.

Size;

S : 34 - 37 number

M : 38 - 41 number

L : 42 - 45 number





Silicone Metatarsal Pad with Toe Separator

Ref: 816

Made of silicone material. This is a metatarsal pillow whose position can be adjusted by placing between toes.

Indications;

- 1-Over-sensitivity of metatarsus head,
- 2-Metatarsalgia,
- 3-Callus developing under the metatarsus,

Size;

Standard Size.



Silicone Metatarsal pad with Loop

Ref: 817

Made of silicone material. Its position is adjusted by putting the silicone ring onto the 2nd toe.

Indications;

- 1-Over-sensitivity of metatarsus head,
- 2-Metatarsalgia,
- 3-Callus developing under the metatarsus,

Size;

Standard Size.



Metatarsal Protective Pad

Ref: 818

This metatarsus pad is made of silicone material; its upper surface is a honeycomb-shaped metatarsal pad with absorbent, comforting and massaging effect. It reduces swelling of the feet and fatigue especially for the people who keep standing for a long time on high heels.

Size;

Standard Size.



Silicone Bunion Shield

Ref: 815

This is a special protective cover made of silicone material which prevent friction on the bunion and protect against trauma. It is used in order to reduce the sensitivity of the bunion.

Size;

Standard Size.



Silicone Toe Ring

Ref: 814

This is a silicone toe cover. It reduces pain and friction on the friction-based bunion.

Indications;

- 1-Toe bunions,
- 2-For reducing the pressure and friction on the toe in deformities such as hammer toe and claw toe.

Size;

S : 34 - 37 number
M : 38 - 41 number
L : 42 - 45 number





Hallux Valgus Shoe

Ref: 821

These shoes help the patient to walk after operatin for hallux valgus deformity. Width is adjustable with the velcro fasteners. Full arch support. The inclination on the sole is effective in reducing the load on the toe.

Size;

Left / Right Available.

S : 35 - 37 number
M : 38 - 39 number
L : 40 - 41 number
XL : 42 - 44 number



Hallux Valgus Shoe (Short)

Ref: 825

Used in the same cases a Ref-821, these shoes provide half-arch-support for hallux valgus.

Size;

Left / Right Available.

S : 35 - 37 number
M : 38 - 39 number
L : 40 - 41 number
XL : 42 - 44 number



Cast Shoe

Ref: 822

When the foot is encased in plaster, the patient can wear these shoes.

Size;

S : 33 - 35 number
M : 36 - 40 number
L : 41 - 45 number



A photograph of a person's legs from the knees down, wearing dark-colored compression stockings. The legs are crossed at the ankles and positioned diagonally across the frame. The background is a smooth gradient from dark brown at the top to a lighter, warm orange-brown at the bottom. The text 'Compression Stockings' is overlaid in the bottom-left corner in a large, white, sans-serif font.

Compression Stockings



Medical Compression Stockings:

- Aurafix medical compression stockings, which have world class compression profile and weaving quality, are a comfortable and quality solution for medical compression stockings.
- Available in CCL 1(18-21 mmHg), CCL2(23-32 mmHg), CCL3(34-46 mmHg) CCL4(49 mmHg and over) compression, open and closed toe.
 - The raw material used is Öko-Tex 100 certificated and minimized allergic reaction of the skin.
 - It carries the moisture of the body outside the stocking, preventing sweating.
 - The elasticity provided by the weaving makes it easy to wear.
 - Due to the strength of the seam, it does not lose its elasticity and durability.
 - With the silicone upper strap which prevents sliding, it can be worn all day comfortably.
 - The closed-toe models have an antibacterial toe weaving, which provides a cleaner and more comfortable use.
 - The open-toe models are soft around the toe, providing easy use and preventing shearing.
 - The heel can carry the load of the body without compromising its elasticity.
 - Medical compression stockings are used in disorders such as chronic venous stasis, lymphedema, post-trombotic syndrome.

	CCL 1 18-21 mmHg Closed Toe	CCL 2 23-32 mmHg Closed Toe	CCL 1 18-21 mmHg Open Toe	CCL 2 23-32 mmHg Open Toe	CCL 3 34-46 mmHg Open Toe	CCL 4 49- mmHg Open Toe
	Product Code No					
A-D Knee-high	AV-101	AV-107	AV-100	AV-106	AV-112	AV-115
A-G Thigh-high	AV-103	AV-109	AV-102	AV-108	AV-113	AV-116
A-T Panty hose	AV-105	AV-111	AV-104	AV-110	AV-114	
A-T/U Pregnancy	AV-117	AV-118				

Thigh-high Medical compression stockings (AV-102, AV-103, AV-108, AV-109, AV-113, AV-116)




Size	Perimeter of the Ankle (cB)	Perimeter of the Lower Calf (cC)	Perimeter of the Upper Calf (cF)
1	18-20 cm	28-33 cm	37-47 cm
2	20-22 cm	30-37 cm	38-52 cm
3	22-24 cm	33-40 cm	43-58 cm
4	24-26 cm	35-43 cm	46-61 cm
5	26-28 cm	37-45 cm	50-66 cm
6	28-30 cm	39-46 cm	53-69 cm
7	30-32 cm	41-48 cm	56-71 cm



Knee-high Medical compression stockings
(AV-100, AV-101, AV-106, AV-107, AV-112, AV-115)

Size	Perimeter of the Ankle (cB)
1	18-20 cm
2	20-22 cm
3	22-24 cm
4	24-26 cm
5	26-28 cm
6	28-30 cm
7	30-32 cm




Compression Stockings



Pantyhose medical compression stockings (AV-104, AV-105, AV-110, AV-111, AV-114)
Pregnancy medical compression stockings (AV-117, AV-118)

Size	Perimeter of the Ankle (cB)	Perimeter of the Lower Calf (cC)	Perimeter of the Upper Calf (cF)
1	18-20 cm	28-33 cm	43-57 cm
2	20-22 cm	30-37 cm	46-60 cm
3	22-24 cm	33-40 cm	49-67 cm
4	24-26 cm	35-43 cm	52-72 cm
5	26-28 cm	37-45 cm	55-76 cm
6	28-30 cm	39-46 cm	60-79 cm
7	30-32 cm	41-48 cm	64-80 cm





Supportive Medical Compression Stockings

This is the medical compression stocking which might be used by those people who work standing, on long journeys, before pregnancy and in order to reduce swelling and pain.

- The elasticity of the weaving provides easy use.
- Patient compliance is high due to the transparency of the weaving.

		A-D Knee-high	A-G Thigh-high	A-T Pantyhose	A-T/U Pregnancy Pantyhose
Product Code No					
140 den	16-18 mmHg	AV-95	AV-96	AV-97	
70 den	12-14 mmHg				AV-98

Pantyhose medical compression (AV-97, AV-98)

Size	Perimeter of the Ankle (cB)	Perimeter of the Lower Calf (cC)	Perimeter of the Upper Calf (cF)
S/1	19-21 cm	28-34 cm	42-57 cm
M/2	22-24 cm	32-38 cm	48-64 cm
L/3	25-27 cm	36-42 cm	54-71 cm
XL/4	28-30 cm	40-46 cm	60-78 cm
XXL/5	31-33 cm	44-50 cm	66-85 cm

Knee-high medical
compression stockings (AV-95)

Size	Perimeter of the Ankle (cB)
S/1	19-21 cm
M/2	22-24 cm
L/3	25-27 cm
XL/4	28-30 cm
XXL/5	31-33 cm

Thigh-high medical
compression stockings (AV-96)

Size	Perimeter of the Ankle (cB)	Perimeter of the Lower Calf (cC)	Perimeter of the Upper Calf (cF)
S/1	19-21 cm	28-34 cm	42-57 cm
M/2	22-24 cm	32-38 cm	48-64 cm
L/3	25-27 cm	36-42 cm	54-71 cm
XL/4	28-30 cm	40-46 cm	60-78 cm
XXL/5	31-33 cm	44-50 cm	66-85 cm



A-D Knee-high	A-G Thigh-high
Product Code No	
AV-99AD	AV-99AG
	AV-99T

Knee-high Anti-Embolism
Stockings (AV-99AD)

Size	Perimeter of the Ankle (cB)
XS	18-20 cm
S	20-23 cm
M	23-26 cm
L	26-29 cm
XL	29-32 cm

Anti-Embolism Stockings


These stockings are used by post-operation inpatients, by women before and after giving birth, by bedridden patients in order to prevent clot forming.

- Comes in 3 models (knee-high, thigh-high, thigh-high with waist strap).
- They have an extra softened section for reducing pressure sore in the heel.
- The observation hole in the metatarsal area provides ease of control.
- The anti-bacterial weaving in the toe section reduces the risk of fungus.
- It has silicone upper straps to prevent sliding.
- The non-latex special weaving assures patient compliance.



Thigh-high Anti-Embolism Stockings
(AV-99AD AV-99T)

Size	Perimeter of the Ankle (cB)	Perimeter of the Lower Calf (cC)	Perimeter of the Upper Calf (cF)
XS	18-20 cm	27-34 cm	37-47 cm
S	20-23 cm	31-37 cm	41-53 cm
M	23-26 cm	35-40 cm	47-58 cm
L	26-29 cm	38-44 cm	52-62 cm
XL	29-32 cm	41-48 cm	58-69 cm





Arm Sleeves

It is used for edama in arms . Comes in 4 different models : CCL1 (18-21 mmHg), CCL2 (23-32mmHg) ve CCL3 (34-46 mmHg) which are compressed and these cover the thumb starting from the wrist.

-It preserves the body temperature due to its high air permeability

-The elastic and soft weaving provides ease of use.




Arm Sleeves (From wrist to shoulder)
(AV-130, AV-131, AV-132, AV-150, AV-151, AV-152)

Size	Perimeter of the Wrist (cC)	Perimeter of the Elbow (cE)	Perimeter of Above Elbow (cG)
1	12-13 cm	17-22 cm	19-25 cm
2	13-15 cm	20-26 cm	22-29 cm
3	15-17cm	23-29 cm	25-32 cm
4	17-19 cm	26-32 cm	29-36 cm
5	19-21 cm	29-35 cm	33-40 cm
6	21-23 cm	32-38 cm	37-44 cm
7	23-25 cm	35-41 cm	41-48 cm






	CCL 1 18-21 mm Hg	CCL 2 23-32 mmHg	CCL 3 34-46 mmHg
From wrist to shoulder (with silicone bandage)	AV-130	AV-131	AV-132
From thumb to shoulder (with silicone bandage)	AV-133	AV-134	AV-135
From wrist to shoulder (with sling)	AV-150	AV-151	AV-152
From thumb to shoulder (with sling)	AV-153	AV-154	AV-155



Arm Sleeves (From thumb to shoulder)
(AV-133, AV-134, AV-135, AV-153, AV-154, AV-155)

Size	Perimeter of Above Palm (cA)	Perimeter of the Wrist (cC)	Perimeter of the Elbow (cE)	Perimeter of Above Elbow (cG)
1	15-16 cm	12-13 cm	17-22 cm	19-25 cm
2	16-19 cm	13-15 cm	20-26 cm	22-29 cm
3	18-21 cm	15-17cm	23-29 cm	25-32 cm
4	20-23 cm	17-19 cm	26-32 cm	29-36 cm
5	22-25 cm	19-21 cm	29-35 cm	33-40 cm
6	24-27 cm	21-23 cm	32-38 cm	37-44 cm
7	26-29 cm	23-25 cm	35-41 cm	41-48 cm



Breast Forms and Bra




Silicone Breast Form

Ref: AG-1050

Easy to use and sleek, the breast form bra is manufactured especially for breast forms. Supporting the breast form, it feels comfortable and provides ease of use. The interior is made of bland cotton fabric. Shoulder straps are softened for easier use. The pockets on both sides allow the breast form to be used both left and right sided.

Size;

15 Sizes (60-130)



Breast Form Bra


Ref: AG-1060

Made of high-quality silicone. It is used in order to capture a natural look after mastectomy, to reduce the pain caused by muscle instability, to protect the operation scar from external impacts and crush.



Size;

15 Boy (60-130)





Rehabilitation Products



Silicone Exercise Ball

Ref: REH-01

These are oval exercise balls made of silicone, used for exercise in rehabilitation. There are 3 hardness levels: Red (hardness), blue (medium), green (soft).



Adjustable Weight

Ref: REH-14

This weight is made of neoprene fabric and consists of pockets to enable the weights to be put on and taken off. There are 1-2-3-4-5 kg options. For 1 and 2 kg options, 250 gr plates are used, while for 3-4-5 kg, 500 gr plates are used. It is used for reinforcing the muscle groups in rehabilitation and for sportsmen.



Silicone Hand Exercise Tool

Ref: REH-15

This is a hand exercise tool made of silicone. It is used in order to prevent lymph edema and to exercise hand muscles during rehabilitation. It comes in standard size and has left & right sides.




Pillows and Seat Cushions

In 1960's, NASA developed the visco elastic memory form in order to reduce the load and pressure on the spinal column of astronauts during their travels. This product has been presented to the industry years later. Our company takes this chance and presents you the aurafix collection, which you will spend 8 hours a day with. These products are manufactured and certified in accordance with European standards.

- The pillows produced using the visco elastic memory form are sensitive to body temperature, pressure and ambient temperature.
- The products support the neck and waist cavity fully.
- It helps your spinal column to remain in a proper line while sleeping.
- It helps you sleep comfortably all night.
- It reduces involuntary motion by 80%.
- Minimizing the pressure, it regulates blood flow and eliminates joint pain.
- It helps removing problems such as stiff neck, lumbar pain and arm numbness.
- It is hygienic, non-allergenic and free of carcinogens.



Pillows and Seat Cushions



Pillows and Seat Cushions





ST-01
(155x54x69 cm)



MS-02
(27x19x34 cm)



MS-01
(21x16x33 cm)



**AT UYGUNLUK BEYANI
EC DECLARATION OF CONFORMITY**

İMALATÇI/MANUFACTURER : TSD SAĞLIK ÜRÜNLERİ TİC. LTD. ŞTİ.

ADRES/ADDRESS : Zümrütevler Mah. Refah sok. Pafti apt. 5A-B
Maltepe / İSTANBUL / TÜRKİYE

TEL/PHONE : +90 216 399 30 55
FAXES / FAX : +90 216 399 63 03

SINIFLANDIRMA
CLASSIFICATION : MDD 93/42/EEC Tıbbi Cihaz Direktifi- EK IX, Kurul 1, SINIF I
Medical Device Directive-ANNEX IX, RULE 1, CLASS I

UYGUNLUK DEĞERLENDİRME : MDD 93/42/EEC Tıbbi Cihaz Direktifi-
EK VII-Uygunluk Beyanı
CONFORMITY ASSESSMENT : MDD 93/42/EEC Medical Device Directive
Annex VII- Declaration of Conformity

ÜRÜN ADI / PRODUCT NAME : EK-1/ANNEX-1

aurafix

TÜRKİYE İLİSTİRLERİNDEN İÇERİLERİN YALNIZCA KULLANILIR TIBBİ CİHAZLAR KATEGORİSİNE GİRER. EEC (93/42) EEC (93/42) GENELİ ŞARTLARA
UYGUN OLUŞUNU BEYAN EDERİZ. DESTEKLEYİCİ TİM BELGELERİ DİREKTİVİN TEGİSLENDİRİLMİŞ BULUNMAKTADIR.
WE HEREBY DECLARE THAT THE ABOVE MENTIONED PRODUCTS MEET THE PROVISIONS OF THE COUNCIL DIRECTIVE 93/42/EEC
FOR MEDICAL DEVICES. ALL SUPPORTING DOCUMENTATION IS RETAINED UNDER THE PREMISES OF THE MANUFACTURER.

İLGİLİ STANDARTLAR /
STANDARTS APPLIED : ISO 13485:2003
MDD COUNCIL DIRECTIVE 93/42/EEC
ISO 14971:2009
EN 980:2008
EN 10411:2008
EN 14155:2011
EN ISO 10993-1:2009
EN ISO 10993-5:2009
EN ISO 10993-10:2009

ÜRETİM TARİHİ
Production Date : 02.01.2008
CE İŞARETİ BAŞLANGIÇ
TARİHİ /
START OF CE MARKING : 02.01.2008 / LOT: 02010800

ŞEHİR, DÜZENLEME TARİHİ
PLACE, DATE OF ISSUE : İSTANBUL / 02.01.2008

İMZA / SIGNATURE : TÜRKAY GÜLHAN
GENEL MÜDÜR
GENERAL MANAGER

